

THE HANDI-PILL ORGANIZER


Thank you for choosing the Handi-Pill Organizer. The Handi-Pill Organizer serves as a perfect device to remind you to take your medication and to store your medication in a safe and handy place. Please read this Instruction Manual thoroughly before using the product. Keep this manual in a convenient place for future reference.

PRODUCT OVERVIEW


Compartment



Compartment Case



LCD Window



Time (Front)

Time (Back)

1. LCD Window

2. Button Cover

3. Minute/Backlight Button
(Min/LT) Set Button (SET)

4. Set Button (SET)

5. Hour Button (HR)

7. Strap Eyelet

8. Clip areas

A. Set Pill Time Indicator

B. Sound Bleep Indicator

C. Clock Indicator

D. Backlight Indicator

E. Time

REPLACING / REMOVING THE COMPARTMENTS FROM THE CASE


Attach the compartment to the timer or the compartment to another compartment using the slots on the side of the compartment or the timer.

INSTALLATION OF THE BATTERIES

The timer unit uses two "AAA" alkaline batteries. To install or replace the batteries, please follow these steps:

1. Remove the battery cover.
2. Insert the batteries observing the correct polarity.
3. Replace the battery cover on the unit.
4. The timer unit is now ready. **Press hold the SET button for 3 seconds to activate the timer**

SETTING THE CLOCK

1. Press the **SET** Key five times until the time and  is blinking.
2. Press the **HR** Button until the correct hour is selected.
3. Press the **MIN/LT** Button until the correct minute is selected.
4. Press the **SET** Key once to confirm.
5. Press **MIN** and **HR** at the same time change 12 /24 format display

SETTING DAILY ALARM

1st time alarm (Morning)

1. Pressing the **SET** key once (while clock is displayed) will display the 1st time alarm setting. The 1 symbol will blink.
2. Press the **HR** Button until the correct hour is selected.
3. Press the **MIN/LT** Button until the correct minute is selected.
4. Press the **SET** Key once to confirm.
5. A dot ● will now appear on top of the 1 symbol.

2nd time alarm (Noon)

1. Pressing the **SET** key twice (while clock is displayed) will display the 2nd time alarm setting. The 2 symbol will blink.
2. Press the **HR** Button until the correct hour is selected.
3. Press the **MIN** Button until the correct minute is selected.
4. Press the **SET** Key once to confirm.
5. A dot ● will now appear on top of the 2 symbol.

3rd time alarm (Evening)


1. Pressing the **SET** key 3 times (while clock is displayed) will display the 3rd time alarm setting. The 3 symbol will blink.

2. Press the **HR** Button until the correct hour is selected.
3. Press the **MIN/LT** Button until the correct minute is selected.
4. Press the **SET** Key once to confirm.
5. A dot ● will be display on top of the 3 symbol.

4th time alarm (Before bed)

1. Pressing the **SET** key 4 times (while clock is displayed) will display the 4th time alarm setting. The symbol will blink.
2. Press the **HR** Button until the correct hour is selected.
3. Press the **MIN/LT** Button until the correct minute is selected.
4. Press the **SET** Key once to confirm.
5. A dot ● will be display on top of the 4 symbol.

RESETTING DAILY ALARM

1. Press the **SET** key until the desired time alarm setting is shown on the LCD display and the appropriate timer symbol (1 or 2 or 3 or 4) is blinking.
2. Press Hour and Min keys at the same time until the alarm is reset to  and the dot ● on top of the timer symbol disappears.
3. Press Set key once to confirm.

LCD BACKLIGHT

Press **MIN/LT** key can turn on the back light, it will last 5 seconds

SPECIFICATION

Timer Unit

- Power Source: 2x AAA Battery operated
- LCD Viewing Area: 16 x 22 mm
- Weight: 29g (w/o Battery)
- Dimension: 28L x 20W x 90H mm (w/o 7 days compartments)
- Battery life for timer unit: Approx. one & a half year

Compartments & Case

- Weight: 96g
- Dimension: 20L x 20W x 90H mm

THE HANDI-PILL ORGANIZER

Instruction Manual