

PRODUCT DESCRIPTIONS & VARIABLES

ROLLZ

Rollz Motion

The Rollz Motion is a walker and wheelchair in one. You can easily transform this walker into a wheelchair and back again in seconds. With the Rollz Motion you can go further and longer by having two mobility aids with you. Features:

- Inverted handles in walker “mode” for optimal posture and less strain on the wrist and joints. The handles can also be adjusted for height.
- Large wheels and a shorter wheel base making for easy steering
- Adjustable lumbar support in the wheelchair seat
- Comfortable seat for resting
- Adjustable footrests
- Collapses easily for transport or storage.

VARIABLES:

- Sizes - Regular ; Small
- Colours: Island Blue (Regular and Small); Matt Black (Regular and Small); Dark Purple, or Pebble White

Rollz Performance

Rollz Motion Performance is the mixed-terrain model of the Rollz Motion walker/wheelchair that comes with air filled tyres for enhanced shock absorption. Other add-ons like the extra soft armrests and wedge-shaped seat add to its comfort.

This walker is designed for frequent use on mixed terrain, for people who might need to sit longer in the wheelchair function or those who need additional cushioning due to sensitive hands, wrists, shoulders or back problems.

VARIABLES:

- Sizes - Regular
- Colours: Jungle Green

SPECIFICATIONS FOR MOTION & PERFORMANCE

	ROLLZ MOTION (SMALL)	ROLLZ MOTION (REGULAR)	ROLLZ PERFORMANCE
Rollator weight	11.6kg	11.6kg	11.4kg
Wheelchair attachment weight	3.4kg	3.4kg	3.4kg

Max weight capacity	125kg	125kg	125kg
User height	1.5m-1.85m	1.6m-1.9m	1.6m-1.9m
Seat height	50cm	55cm	55cm
Handle height	82cm-94cm	87cm-98cm	87cm-98cm
Seat size (rollator)	20cmx45cm	20cmx45cm	20cmx45cm
Seat size (wheelchair)	42cmx45cm	42cmx45cm	42cmx45cm
Overall width	67cm	67cm	67cm
Overall length	63cm	63cm	63cm
Width when folded	31cm	31cm	29cm
Front wheel diameter	8 inch	12 inch	8x 1.25 inch
Back wheel diameter	8 inch	12 inch	12x 2.5 inch

Rollz Motion/Performance Cup Holder

The cup holder slides easily on the frame and allows to transport a bottle or cup. It can remain in place when folding the walker. The Rollz Motion cup holder is only intended for use with the Rollz Motion and Performance

Rollz Motion/Performance Travel Cover

Protect your walker during air or bus trips with the sturdy travel cover. When folded, the walker and the wheelchair package are both protected in the compact cover. Includes a wheelchair label that can be easily attached to the travel cover. This accessory fits the Rollz Motion, Rollz Performance and Rollz Motion Rhythm.

Rollz Motion/Performance 3-in-Wheelchair Package Holder

This 3-in-1 holder includes two frame clips and one cane clamp and can hold the wheelchair package, a cane and a bag. The clips can be attached to the frame and provides an alternative place to carry the wheelchair package above the seat, for more walk-in space. The clamp can be attached to the cane and than hanged on the outside of the holders. An extra bag can also be attached on the holders. The cane clamps come with three different elastic attachments to fit most cane sizes. This holder can be used with all Rollz Motion and Performance models.

Rollz Motion/Performance Back Support

This flexible back support is an useful accessory when you go out and use your Rollz Motion or Performance as a walker. You can lean against it while sitting on the rollator, when you

do not have the wheelchair package with you. The back support can be attached to the frame only when the Rollz Motion/Performance is used in rollator position.

Rollz Motion/Performance Seat Belt

The seatbelt can be fixed with Velcro at the back of the wheelchair package and adjust to the desired length. Make sure the wider, foamed part is at the front of your abdomen. When the seatbelt is not in use, you can unclip the buckle and store it at the top of the wheelchair package. You can convert the Rollz Motion into a walker while the belt is attached to the wheelchair package. Works with all the Rollz Motion and Performance models. Note: The seat belt is not a safety belt.

Rollz Motion/Performance Umbrella

The umbrella protects the user from rain and sun and can be attached to the frame so the rollator can still be used with both hands. The umbrella is adjustable so it can be angled in different ways. This accessory is suitable for all the Rollz Motion/Performance models.

Rollz Flex

Rollz Flex is a compact and stable walker with a large shopping bag. Its innovative design enables the walker to stay close to the body so that you can walk with optimum posture. The short wheel base makes it easy to step aside from it without the wheels getting in the way. It also allows you to make tight turns in narrow spaces.

This walker offers great support and steers lightly, even with a full shopping bag which can carry a load of up to 20kg without affecting the stability of the walker. The internal bag inside the carrier enables you to easily lift the contents out after a shopping trip. The softer foam filled tyres of the Rollz Flex 2 absorb vibrations, allowing a more comfortable walk. The beautiful design hides the more traditional walker look.

VARIABLES:

Size: Regular

Colour: White Frame with Grey Bag

Technical Data	REGULAR
Weight	7.4kg
Handle Height	80 – 100 cm
Width between handles	58cm
Seat width	44cm
Seat height	62cm
Width	61cm
Width when folded	33cm
Length	60cm

Rollz Flex Cane Holder

The cane holder consists of two parts: a frame clip and a cane clamp. The frame clip is slid onto the frame of the Rollz Flex and the clamp is attached to the cane with a supplied elastic. You can easily slide your cane into the frame clip and remove it when you want to use it. The cane holder fits different stick sizes thanks to the multi-size elastics that are included in the package.

Rollz Flex Shopping Hooks

This set of two shopping hooks allows to bring along extra bags with light items. The shopping hooks are easy to attach to the handle bar of the Rollz Flex without tools. This way you can keep both hands on the handle bar and let the walker carry the shopping bags for you.

byACRE

byACRE Carbon Ultralight

Perfect for adventures around the world or just around the corner.

- Only 4.8 kg and extremely easy to fold
- Carbon fibre frame
- Ergonomic, reversed handles for optimal posture and reduced pressure on wrists and joints
- Easy height adjustment
- Comes fully assembled

VARIABLES:

Colour: Strawberry Red; Carbon Black or Oyster White

Sizes: Regular or Wide Track

SPECIFICATIONS:

Technical Data	Regular	Wide Track
Weight	4.8kg	4.95kg
Handle Height	79-91cm	79-91cm
Width between handles	39cm	44cm
Seat width	42cm	47cm
Seat height	61cm	61cm
Width	62.5cm	66.5cm
Width when folded	25cm	25cm
Length	66cm	66cm

byACRE Carbon Overland

Carbon Overland walker is a sturdy and versatile mixed-terrain walker boasting the looks of the latest SUV.

> Lightest mixed-terrain rollator at 6.7 kg for easy lifting in your vehicle

- > Foldable, fully assembled, to make it hassle free to start exploring
- > Extremely stable carbon fibre frame – 5 x stronger than steel
- > Ultra-durable pneumatic rubber tires for better suspension
- > Beach and snow-friendly with rigid front forks and sealed bearings
- > Aircraft grade aluminium in gunmetal colour, for strength and lightness

VARIABLES

Colour: Green or Brown

Size: One size

Weight	6.7 kg
Handle height	82-95 cm
Seat width	47 cm
Seat height	61 cm
Width	66 cm
Width, folded	25 cm

byACRE Organiser Bag

Need a small bag for your walker? The organiser is designed to hold your personal belongings such as keys, wallet, phone. Compatible with both the Ultralight and the Overland (wide track for Overland).

Variations: Regular and Wide Track

byACRE Back Support

To better facilitate a rest while out on a trip, a backrest that supports your lower back while sitting on the seat of your rollator. Compatible with both the Ultralight and the Overland (wide track for Overland).

Variations: Regular and Wide Track

byACRE Weekender Bag

Going on a weekend trip or small adventure? The weekend bag offers a convenient way to carry your necessities. Compatible with both the Ultralight and the Overland (wide track for Overland).

byACRE Cane Holder

To hold your cane in place while walking with your Carbon Ultralight rollator we have designed this smart and easy-to-use cane holder. Compatible with both the Ultralight and the Overland (wide track for Overland).

byACRE Ultralight Carry Bag

- Wear-resistant, padded and sturdy material with an elegant and soft lining inside for optimal protection
- Carry-handle and shoulder strap for easy transportation in airports, busses, trains or cruise ships
- Integrated protection pad so the wheels and axles are protected during transport
- Perfect custom fit for the best protection possible

Only compatible with the byACRE Carbon Ultralight.