

# KIT EAR CUSHION

References	Designations
NEG1103	KIT EAR CUSHION 30 x 30 cm
NEG1103HLEN	LENZING COVER FOR EAR CUSHION 30 x 30 cm - SET OF 2



- Supports the patient's head
- Relieves the ears and prevent bedsores

## COMPOSITION :

- Cushion's protective cover: bi-elastic fabric, with a 54% coating of polyurethane, 46% polyamide
- Cushion's filling: 100% polystyrene micro-beads
- Bi-elastic PU cover: coated with 54% of polyurethane, 46% polyamide
- Lenzing FR® cover: 70% lyocell, 30% viscose coating with 100% of polyurethane

## OVERALL RECOMMENDATIONS

- Control regularly the overall state of the cover (wear and tear, etc...) and do not use the cushion anymore if the cover is faulty.
- In order not to damage the cover, the use of solvent or chlorine-based cleaning products are not recommended. Instead, we strongly encourage to mix chlorine-based cleaning products with water and rinse several times.
- The sterilization of the cover is not recommended and could damage the product.

## MAINTENANCE :

Removable insert: clean with wipes

Do not machine wash (the protective cover could crack and it may damage the washing machine)

Cushion :



Bi-elastic PU cover & Lenzing FR® cover:

