




# BEDROOM SAFETY & SUPPORT

The Aspire bedroom safety & support range has been designed with user support and safety in mind. Bed poles, cradles and rails ensure carers can provide safety and support for residents without compromising comfort in their bedroom environment.





Drop Down Bed Rail BEA013700



Bed Pole With Adjustable Crook Handle BEA014150



Bed Pole With Adjustable Crook Handle BEA013100



Fixed Bed Rail BEA013600



Clamp On Bed Pole BEA009300 & BEA009310



Adjustable Height Bed Cradle



Fixed Height Bed Cradle





# **aspire** BEDROOM SAFETY & SUPPORT



#### FREE STANDING SELF HELP POLE

#### BEA013100

- Free standing support handle provides assistance for the user to reposition or transfer in and out of bed
- · Handle height is easily adjusted to suit the user
- · Sits behind or beside the bed





# BED POLE WITH ADJUSTABLE CROOK HANDLE

#### BEA014150

- $\cdot$   $\,$  Support pole designed to assist the user to balance while getting into and out of bed
- · Adjustable crook handle can be located in multiple positions
- · Base is positioned under the mattress



#### **CLAMP ON BED POLE**

#### BEA009300 (RIGHT) & BEA009310 (LEFT)

- Support pole clamps on to most bed frames to provide assistance with repositioning or transferring in and out of bed
- · Folding crook handle allows the pole to swing down and out of the way as required
- · Clamps onto the side of most bed frames



#### **FIXED BED RAIL**

# BEA013600

- Support and safety rail provides repositioning and transfer support and assists to prevent falls from bed
- · Base is positioned under the mattress



# **DROP DOWN BED RAIL**

#### BEA013700

- · Support and safety rail when raised is effective at preventing falls from bed
- · May be lowered down and out of the way for patient transfers or repositioning
- · Base is positioned under the mattress



## **FIXED HEIGHT BED CRADLE**

### BEA014605

- Keeps bed linen off sensitive areas of the patient by creating clearance between linen and mattress
- · Base is positioned under the mattress



# ADJUSTABLE HEIGHT BED CRADLE

# BEA014610

- Keeps bed linen off sensitive areas of the patient by creating clearance between linen and mattress
- · Adjustable height suitable for a range of mattress heights
- · Base is positioned under the mattress