

BACK CUSHION XL

References	Designations
NEG1301	Back cushion XL
NEG1436	Loop cotton cover
NEG1430	Bi-elastic PU cover
NEG1301HLEN	Lenzing FR® cover
Dimensions : 260 x 60 cm	



This cushion supports the patient in sitting position or in lateral decubitus position in bed.
Depending on types of armchairs, it can be used as a wedge


COMPOSITION :

- Cushion's protective cover: bi-elastic fabric with a 54% of polyurethane coating, 46% polyamide
- Cushion's filling: 100% polystyrene micro-beads
- Loop cotton cover: 90% cotton, 10% polyester
- Bi-elastic PU cover: 54% of polyurethane coating, 46% polyamide
- Lenzing FR® cover: 70% lyocell, 30% viscose coated with 100% of polyurethane

OVERALL RECOMMENDATIONS :

- Control regularly the overall state of the cover (wear and tear, etc...) and do not use the cushion anymore if the cover is faulty.
- In order not to damage the cover, the use of solvent or chlorine-based cleaning products are not recommended. However, if it is necessary to use chlorine-based cleaning products, we strongly encourage to mix it with water and rinse several times.
- The sterilisation of the cover is not recommended and could damage the product

MAINTENANCE :

Cushion : 



Loop cotton cover : 



Bi-elastic PU cover & Lenzing FR® cover : 