

# DECUBITUS LATERAL CUSHION XL

References	Designations
NEG0610	Decubitus Lateral cushion XL
NEG0610D	Decubitus Lateral cushion XL with removable cover
NEG1413	Interlock cover
NEG0631	Loop cotton cover
NEG0624	Bi-elastic PU cover
NEG0610HLEN	Lenzing FR <sup>®</sup> cover
<b>Dimensions : 200 x 35 cm</b>	



This cushion provides a great support for the whole body, from head to feet. When in decubitus lateral position, it relieves the pressure effectively and prevent frictions around the arms, knees, and calves. In any lateral position, it is a great positioning cushion.

## COMPOSITION :

- Cushion's protective cover: bi-elastic fabric with 54% of polyurethane coating, 46% polyamide
- Cushion's filling: 100% polystyrene micro-beads
- Removable insert: 100% cotton mesh Jersey
- Interlock cover : 100% cotton
- Loop cotton cover: 90% cotton, 10% polyester
- Bi-elastic PU cover: 54% of polyurethane coating, 46% polyamide
- Lenzing FR<sup>®</sup> cover: 70% lyocell, 30% viscose coating with 100% of polyurethane

## OVERALL RECOMMENDATIONS

- control regularly the overall state of the cover (wear and tear, etc...) and do not use the cushion anymore if the cover is faulty.
- In order not to damage the cover, the use of solvent or chlorine-based cleaning products are not recommended. However, if it is necessary to use chlorine-based cleaning products, we strongly encourage to mix it with water and rinse several times.
- The sterilization of the cover is not recommended and could damage the product.

## MAINTENANCE :

Removable insert: cleaning by using wipes

Do not machine wash (the protective cover could crack and it may damage the washing machine)

Cushion :

Interlock cover :

Loop cotton cover:

Bi-elastic PU cover & Lenzing FR<sup>®</sup> cover: