


Instruction Manual of Pill Box reminder

HH !'

Features

- Displays current time – hour and minute
- Power Saving Alert (PSA), continuously reminding three times every minute forever unless MODE/CONFIRM key is depressed.
- Up to four reminders a day
- Use 1 pc of LR44 battery



Keys

1. MODE – Used to enter each reminder mode, to return to real time, or to confirm dose taken, and stop the reminder alert
2. SET – To enter programming mode of real time or each reminder alarm time
3. ADJ – To adjust real time, and each reminder alarm time, or to enable and disable reminder.



How to Set Real Time

1. At real time mode press "SET" key to set the hour, the reading of "hour" starts to flash. Use "ADJ" key to adjust the hour.
2. Press "SET" key again to set the minute, the reading of minute starts to flash. Use "ADJ" to adjust the minute.
3. Press "MODE" again to finish time setting and return to real time.

How to Set Reminder Time

At real time mode, press the 'MODE' button once to scroll to the 1st Reminder Alarm setting. The Reminder Mode symbol displays solid with 1st reminder symbol flashing, indicating it is not yet enabled, as shown in Fig a. To adjust the reminder time, use the same method you used when you set the real time, as in the above paragraph. Once set, make sure the "Hour" and "Minute" reading are displayed solid. If the 1st reminder symbol is flashing (refer Fig. a). This means the 1st reminder is not yet enabled. Press "ADJ" to enable the reminder. The symbol becomes solid, indicating it is enabled. If required, press "ADJ" again to disable the reminder. The symbol starts to flash again. Use "MODE" again to scroll to the 2nd, 3rd, or 4th reminder, and apply the same method to set each reminder.

How to Know the Reminder is Turned On

At real time, if the reminder symbol appears for a certain reminder, as shown in Fig a., then that reminder is turned on. If the symbol does not appear for a certain reminder, then that reminder is turned off.

How to Enable or Disable Each Reminder

At real time mode, repeatedly press "MODE" to scroll to the reminder, which you are going to enable. If the reminder symbol is flashing then it is not enabled. Press "ADJ" to enable it. The related symbol should become solid.

How to Confirm the Reminding Alert

The symbols of all enabled reminders should display solid at real time mode. When the reminder time arrives, its related symbol starts to flash and the reminding alert goes off for three times, then repeats every minute. You don't need to open the cover, just press the red button to confirm dose is taken, and stop the alert. Otherwise the alert will never stop. After red button is pressed the symbol becomes solid again. It will repeat again the next day at the same time. Be sure to use the same method to confirm the reminder and stop the alert. If your need for a reminder becomes unnecessary, then follow the instruction above to disable the reminder.

How to Replace Battery

One piece of LR44 is used. Please follow the steps to replace battery:

1. Use screwdriver to remove four screws on the backside.
2. Use fingers to remove the front case carefully. Turn it upside down.
3. Use screwdriver again to remove two screws "A", and then take out battery stopper "B".
4. Remove battery "C" as shown. Discard the old battery.
5. Insert a fresh battery (positive side up) into battery holder "D".
6. Put back battery stopper "B", and screw back the two screws "A".
7. Install back the front case to main case properly. Make sure the buzzer "F" seats on its holder "E" properly.
8. Use driver again to screw back the four screws.

