

# Mattress Selector




Involuntary movement during sleep  
[approx. every 15 mins]


PATIENT MOBILITY  
TISSUE DAMAGE  
PRODUCT OPTIONS

Able to turn and move independently or with minimal assistance in bed.

No damage



Prevention



**Flex**  
SWL 213kg / Clinical SWL 160kg

Restricted ability  
Some assisted repositioning required.

Grade 1-2 pressure damage or Grade 4 heel pressure damage



Level 1-2 Pressure Wound



**Aria+** SWL 230kg




**Fika** SWL 159kg

Does not move during sleep & / or spends all in bed. Minimal ability to move in bed, Assisted repositioning required on most occasions..


Grade 4 pressure damage or grade 3 pressure damage in a patient with a high/low body mass index



Level 3-4 Pressure Wound



**Aria+ with Pump** SWL 250kg



**Promatt** SWL 227kg



**Alto** SWL 175kg