

IC30 ACTIVEX™ MATTRESS



The IC30 is our all in one mattress. Incorporating pressure care, bed mobility and most importantly a very comfortable sleeping surface. Highly effective in the treatment and prevention of pressure sores.



RISK CATEGORY:

HIGH RISK

COMFORT CATEGORY:

MEDIUM FEEL



Aged Care




Every Day Care



Disability



Facilities



Cooling
Comfort
Care


SLEEP QUALITY

Quality sleep is an important factor to wellbeing and healing. The IC30 mattress creates an ideal balance for medical and pressure care needs, sleep quality and longevity.



TRANSFERS

Transferring to and from the IC30 mattress is assisted with the XP Support mattress edge and anti-hear BlueTec mattress surface cover.



PRESSURE CARE

The IC30 Mattress is designed to prevent and reduce the risk of pressure injuries. The ActiveX™ sleep surface provides partial immersion with envelopment, creating an excellent average peak pressure, while maintaining bed mobility and lateral movement. The mattress microclimate is regulated by using BlueTec, gel infusion and charcoal.


FEATURES




Pressure Relieving




ZeroTec covers available



ActiveX™ Gel temperature regulating



Anti-shear materials



Moisture wicking & breathable fabric.




ActiveX™ Charcoal reduces odors and is microbial.



Multi layer immersion



Reduces Rollout Risk




Contour-Latch™ attaches to most profiling beds.



20+ waterlow rating


BED CODES & SIZES




SINGLE
Code: IC30S
Mattress Size: 1900 x 900




LONG SINGLE
Code: IC30LS
Mattress Size: 2030 x 900



KING SINGLE
Code: IC30KS
Mattress Size: 2030 x 1070




DOUBLE
Code: IC30D
Mattress Size: 1900 x 1350



LONG DOUBLE
Code: IC30LD
Mattress Size: 2030 x 1350



QUEEN
Code: IC30Q
Mattress Size: 2030 x 1520



SPLIT QUEEN
Code: IC30SQ
Mattress Size: 2030 x 750 (per side)

SPECIFICATIONS

ActiveX™ pressure relieving material
Vant ventilated sides and base fabric
Bluetec temperature regulating cover
Option Zerotec waterproof cover (see pg. 22)
Flip free design
Zip off machine washable covers
30kg - 200kg user weight
Standard cover MVTR: NA.
Medical Cover MVTR (g/24h·m²): 1635.0 Average
Medical Cover BS7175: Crib 5
Waterlow ≤20

Feel Rating	Medium
Mattress Height	200mm
Risk Category	High Risk
User Weight	20kg - 228kg
Warranty	10 Year Warranty

PRESSURE CARE

Risk Category High

This rating is based on our average peak pressure of <19mmHg. For some high risk clients Icare Zerotec covers are recommended.

Icare risk category ratings are based on average peak pressure, Therapist recommendation, case studies, industry comparison and the Waterlow scale. For some high risk clients Icare Zerotec covers are recommended. The best way to assess a client on a pressure reducing surface is to monitor pressure points with visual checks over a period of 5 days.

CARE & CLEANING

To launder the cover in a washing machine, set the main wash to a minimum of 65°C for no less than 10 minutes. The cover can be tumble dried on low heat and the temperature must not exceed 60°C. The foam can be autoclaved at 134°C.

For Sales: **1300 848 077**

backtosleep.com.au